

STRONG & PRECISE


LAPAROSCOPIC INSTRUMENTS

Superior Shaft Strength

37% Stronger than other instruments
in the market

Lowest shaft deformity within long surgery
in comparison to other instruments
in the market

Highest time and load durability



Developing Products in the Field of Laparoscopic Surgical
Instruments and Continually Plan and Develop Products
that Actively Represent the Ideas of the Excellent Medical
Personnel and Voices of the Field.



ACH Medical
Achieving the Future

ACH Medical Co., Ltd.

28, Namyang-ro 930beon-gil, Namyang-eup,
Hwaseong-si, Gyeonggi-do, Rep. of Korea, 18255
TEL +82-31-355-7110 FAX +82-31-355-6759

Copyright ACH Medical.,Ltd. All Rights Reserved.

NL-LI-2004-REV03



LAPAROSCOPIC HAND INSTRUMENTS



Optimized Force Transmission for Grasping and Dissecting
Enables maximum force transmission during closure of jaws, ensuring strong and precise control



Selection of the Most Commonly Used Jaws for Single Use
Scissors, Dissector, Fine dissector, Forceps, Long Forceps, Corce-olmi, Clinch, Crocodile forceps, Babcock forceps



Classification Tips by Different Color 3 Jaw with Each Color
Color coding allows clear identification of the scissors grasping Forceps, dissector



Smoother Rotation of Knob
Comfortable single-handed rotation knob



The Longest Laparoscopic Instruments World-wide
Available in 33cm, 36cm, 45cm and 47cm shaft lengths

SCISSORS | Ø5mm



SCISSORS

METZENBAUM, Multipurpose Curved Scissors

	33cm	36cm	45cm	47cm
Lock-free	LHSC 0533FN	LHSC 0536FN	LHSC 0545FN	LHSC 0547FN
Ratchet				



SCISSORS, INSULATED TO TIPS

METZENBAUM, Multipurpose Curved Scissors

	33cm	36cm	45cm	47cm
Lock-free	LHSC 0533F	LHSC 0536F	LHSC 0545F	LHSC 0547F
Ratchet				

DISSECTOR | Ø5mm



DISSECTOR

MARYLAND, Multipurpose Curved Dissection Forceps

	33cm	36cm	45cm	47cm
Lock-free	LHDB 0533F	LHDB 0536F	LHDB 0545F	LHDB 0547F
Ratchet	LHDB 0533R	LHDB 0536R	LHDB 0545R	LHDB 0547R



FINE DISSECTOR

MARYLAND,
Needle Holder Assistance and Delicate Dissection Forceps

	33cm	36cm	45cm	47cm
Lock-free	LHDS 0533F	LHDS 0536F	LHDS 0545F	LHDS 0547F
Ratchet	LHDS 0533R	LHDS 0536R	LHDS 0545R	LHDS 0547R

GRASPING FORCEPS | Ø5mm



FORCEPS

Multipurpose Grasping Forceps, Fenestrated

	33cm	36cm	45cm	47cm
Lock-free	LHFC 0533F	LHFC 0536F	LHFC 0545F	LHFC 0547F
Ratchet	LHFC 0533R	LHFC 0536R	LHFC 0545R	LHFC 0547R



LONG FORCEPS

Multipurpose Grasping Forceps, Fenestrated

	33cm	36cm	45cm	47cm
Lock-free	LHLF 0533F	LHLF 0536F	LHLF 0545F	LHLF 0547F
Ratchet	LHLF 0533R	LHLF 0536R	LHLF 0545R	LHLF 0547R



CROCE-OLMI FORCEPS

Manipulation of Tissue Grasping Forceps

	33cm	36cm	45cm	47cm
Lock-free	LHCO 0533F	LHCO 0536F	LHCO 0545F	LHCO 0547F
Ratchet	LHCO 0533R	LHCO 0536R	LHCO 0545R	LHCO 0547R



CLINCH FORCEPS

Atraumatic Grasping Forceps


	33cm	36cm	45cm	47cm
Lock-free	LHCL 0533F	LHCL 0536F	LHCL 0545F	LHCL 0547F
Ratchet	LHCL 0533R	LHCL 0536R	LHCL 0545R	LHCL 0547R



CROCODILE FORCEPS

Strong Fixation of Tissue Grasping Forceps

	33cm	36cm	45cm	47cm
Lock-free	LHCC 0533F	LHCC 0536F	LHCC 0545F	LHCC 0547F
Ratchet	LHCC 0533R	LHCC 0536R	LHCC 0545R	LHCC 0547R



BABCOCK FORCEPS

Manipulation of Large Fragile Tissue Grasping Forceps

	33cm	36cm	45cm	47cm
Lock-free	LHBC 0533F	LHBC 0536F	LHBC 0545F	LHBC 0547F
Ratchet	LHBC 0533R	LHBC 0536R	LHBC 0545R	LHBC 0547R