

COMFORT & EASY CONTROL

LAPAROSCOPIC SUCTION & IRRIGATION SYSTEM



LAPAROSCOPIC SUCTION & IRRIGATION SYSTEM



Highest Suction Performance

Shaft with continuous inner sheath prevents clogging



Simplified Structure

Simplified internal mechanism minimizes the chance of leaks



Variety of Disposable Probe Length and Tips

Variety of disposable probe lengths(29, 34, 36, 45cm) and wide range of tip offerings



Soft Co-joined Tube

Low torque tubing with easy separation



Ergonomic Design

Improved grip handle for optimized ergonomics

PRODUCT CODE



Electrode Option

BASIC

Shaft Diameter Ø5mm

	29cm	34cm	36cm	45cm
Tubing	LSBHI0529T	LSBHI0534T	LSBHI0536T	LSBHI0545T
Connector	LSBHI0529B	LSBHI0534B	LSBHI0536B	LSBHI0545B



Electrode Option

L-HOOK

Shaft Diameter Ø5mm

	29cm	34cm	36cm	45cm
Tubing	LSLHI0529T	LSLHI0534T	LSLHI0536T	LSLHI0545T
Connector	LSLHI0529B	LSLHI0534B	LSLHI0536B	LSLHI0545B



Electrode Option

J-HOOK

Shaft Diameter Ø5mm

	29cm	34cm	36cm	45cm
Tubing	LSJHI0529T	LSJHI0534T	LSJHI0536T	LSJHI0545T
Connector	LSJHI0529B	LSJHI0534B	LSJHI0536B	LSJHI0545B



Electrode Option

SPATULA-CURVED

Shaft Diameter Ø5mm

	29cm	34cm	36cm	45cm
Tubing	LSSCI0529T	LSSCI0534T	LSSCI0536T	LSSCI0545T
Connector	LSSCI0529B	LSSCI0534B	LSSCI0536B	LSSCI0545B



Electrode Option

SPATULA-STRAIGHT

Shaft Diameter Ø5mm

	29cm	34cm	36cm	45cm
Tubing	LSSSI0529T	LSSSI0534T	LSSSI0536T	LSSSI0545T
Connector	LSSSI0529B	LSSSI0534B	LSSSI0536B	LSSSI0545B